

# Liteway 8



**Transportable scooters**

**Your new route to independence**

# Liteway 8

## Liteway 8 Specification:

Type Class	C / 3
Overall length	120cm / 47.2"
Overall width	59cm / 23"
Floor to seat height	63cm (min)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	3 sets of 12V 14Ah x 2
Maximum speed	12.4km/h / 7.8mph
Ground clearance	8cm / 3.2"
Maximum kerb height	5cm / 2"
Range*	up to 38km / 23.8 miles
The maximum safe slope	8°
Turn-around width	1.40m / 4'6"
Colour	Blue / Graphite

\*Range on full charge and flat ground.

## Main Features:

- Car transportable
- Rear coil spring suspension
- Flip-up width adjustable armrests
- Handy storage compartment
- Tiller adjustment and can be folded down for ease of transportation
- Swivel seat
- Single heaviest component (rear chassis) 18.1kgs / 40lbs

## Optional Extras:

- Optional accessory a front basket
- Optional lighter seat
- Optional upgrade of an off-board charging station



## Liteway 8 Profile:

At last, a Class 3, transportable, aluminium scooter that has been especially designed to go the extra mile. The Liteway 8 provides the benefits of an 8 mph scooter in terms of comfort and range of travel, yet offering the convenience of fitting easily into the boot of most family-sized cars.



Lightweight Components



Easy to Dismantle



User-friendly Dashboard



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E. - v 2 Final-August 2014