



Rascal 388 S & 388 Deluxe



*Photograph shows 388 S model in Silver & 388 Deluxe model in Red.



Mid-size scooters

Your new route to independence

Rascal 388 S & 388 Deluxe

Rascal 388 S and 388 Deluxe Specification:

Type Class	B / 2
Overall length	125cm / 49"
Overall width	54cm / 21"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 33Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	10cm / 4"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	2.87m / 9'5"
Colour 388 S	Red / Silver
Colour 388 Deluxe	Red

*Range on full charge and flat ground.

Main Features 388 S:

- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis

Main Features 388 Deluxe:

- All of the above features of the 388 S
- Lift-up tiller head for easy access
- Upgraded captains seat
- Smooth controls

Optional Model:

- The Rascal 388 XL, offers all the great features of the 388 Deluxe, plus a maximum speed of 6 mph and a ground clearance of 13cm



Rascal 388 S and 388 Deluxe Profile:

The Rascal 388 S and 388 Deluxe are both well established, great value for money and suit your everyday mobility needs. At 21" wide, it is the narrowest width of any scooter in its class it is equally at home inside a busy supermarket as it is on the pavement. With a comfortable high back seat with headrest, smooth controls and adjustable tiller, these high performance, reliable scooters are your ideal mobility partner.

Floor to seat height on the 388 S model is 56-64cm (min-max), seat width 45.7cm / 18" seat depth 40.6cm / 16".

Floor to seat height on the 388 Deluxe model is 57-65cm (min-max), seat width 45.7cm / 18" seat depth 45.7cm / 18"



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 3 Final-August 2014

