


Rascal Vantage X



Delta Tiller



Waterproof Panel and Twin Mirrors



Maximum Speed



Rascal Vantage X

Rascal Vantage X Specification:

Type Class	C / 3
Overall length	120cm / 47.2"
Overall width	58.5cm / 23"
Floor to seat height	53-58cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	149.7kgs / 23.5 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	6.5cm / 2.6"
Maximum kerb height	7.6cm / 3"
Range*	up to 42km / 26 miles
The maximum safe slope	10°
Turn-around width	1.40m / 4'6"
Colour	Black / Red

*Range on full charge and flat ground.

Main Features:

- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Infinite adjustable delta tiller position
- Modern delta handlebars
- Waterproof panel
- Emergency hydraulic braking system



Rascal Vantage X Profile:

The Rascal Vantage X is a compact, value for money 6 mph road scooter with extremely smooth user-friendly controls. With a comfortable seat, infinitely adjustable delta-bar tiller steering and waterproof switches this highly manoeuvrable scooter features LED lighting and can adequately accommodate users up to 23.5 stone.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 4 Final-August 2014