

Rascal 850



Performance scooters

Your new route to independence

Rascal 850

Rascal 850 Specification:

Type Class	C / 3
Overall length	127cm / 50"
Overall width	65cm / 26"
Floor to seat height	62cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	15cm / 6"
Maximum kerb height	10cm / 4"
Range*	up to 38km / 24 miles
The maximum safe slope	9°
Turn-around width	3.12m / 10'2"
Colour	Black / Red

*Range on full charge and flat ground.


Main Features:

- Five-position flip-up tiller
- Sliding luxurious leather seat
- Flip-up width adjustable armrests
- All round coil suspension
- Height adjustable headrest
- Backrest angle alteration
- Sturdy front bumper



Rascal 850 Profile:

Whether driving at 4 mph on pavements or 8 mph on roads, the Rascal 850 brings together compact design with Class 3 specification. Available in two colour options and with 50Ah batteries as standard you can travel in confidence and style. A fully adjustable sculpted sliding seat with added lumbar support combined with full suspension ensures a smooth ride, whatever the terrain.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 3 Final-August 2014