

RASCAL

**Performance scooters** 

Your new route to independence

# **Rascal Ventura**

## **Rascal Ventura Specification:**

TO	
Type Class	C/3
Overall length	131cm / 51.6"
Overall width	68.4cm / 26.9"
Floor to seat height	63-68cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Max carrying capacity - batt pk (50Ah) <sup>†</sup>	159kgs / 25 stone
Battery voltage and capacity**	12V / 50Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range*- standard battery pack (50Ah)	up to 37km / 23 miles
Range*- upgraded battery pack (75/80Ah)	up to 59km / 36 miles
The maximum safe slope	10°
Turn-around width	1.85m / 6'
Colour	Black / Silver
*Range on full charge and flat ground. **Upgrade to 12V / 75/80Ah x 2 <sup>†</sup> Max carrying capacity - batt pk (75/80Ah) 133.8kgs / 21 stone	

### **Rascal Ventura Profile:**

#### Main Features:

- Front coil and rear mono-shock suspension for a superior ride
- Gas strut tiller with easy access battery charging port
- Modern delta handlebars
- Modern LCD panel and twin mirrors
- Safety turning sensor, slows down scooter when turning at full speed
- Auxiliary USB power outlet for charging mobile phones etc
- Emergency braking system
- Modern LED lighting system

## **Optional Extras:**



- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- Standard rain canopy upgrade to a polycarbonate front shield for an additional cost

The Rascal Ventura is a robust, comfortable 8 mph scooter that will take you wherever you wish to go. With its modern practical design and state of the art styling the Ventura is a real head turner. Couple this with its delta handlebar, tiller adjustment, LED lighting, LCD instrument panel and a waterproof USB power outlet the Ventura has all the makings of a contemporary classic.

