

Technical Specification - Ultralite 765

Designation	Model 765	Model 760
Dimensions (L x W x H) ¹	38" x 22" x 38" (96 x 57 x 96 cm)	38" x 22" x 38" (96 x 57 x 96 cm)
Weight Capacity ²	200 lb. (91 Kg)	200 lb. (91 Kg)
Maximum Speed	3.8 mph (6 km/h)	3.8 mph (6 km/h)
Safe Incline	6° (10%)	6° (10%)
Max. Incline ³	10° (17%)	10° (17%)
Range ⁴	9 miles (14 km)	9 miles (14 km)
Turning Radius	30" (77 cm)	30" (77 cm)
Total Weight	99 lb. (45 Kg)	88 lb. (40 Kg)
Total Weight w/out batteries	77 lb. (35 Kg)	66 lb. (30 Kg)
Seat	Detachable & Folding	Detachable & Folding
Battery Capacity	12V/12Ah (2)	12V/12Ah (2)
Battery Pack Weight	22 lb. (10 Kg)	22 lb. (10 Kg)
Charger	Off Board (2A)	Off Board (2A)
Motor Output Power	0.45 hp / 336W	0.45 hp / 336W
Front Tyre Size	6" (15.34 cm)	6" (15.34 cm)
Rear Tyre Size	8" (20 cm)	8" (20 cm)

Notes:

1. Including anti-tip casters
2. On level terrain
3. Based on rider weight of 165 lb.
4. Range varies with rider's weight, terrain, battery charge, temperature and style of driving