

Rascal P320 Compact



Up to 20km
12.5 Mile
Range*

A black icon of a winding road with white dashed lines, representing the range of the powerchair.

Great Turn-around Width
97cm / 3'2"

A white arrow forming a U-turn, with the text 'Great Turn-around Width' written along its curve.

Max. Carrying
Capacity
136kgs / 21 Stone

A white icon of a speedometer with a needle pointing to the right.

Easy to
Dismantle

A black puzzle piece and a white puzzle piece, one of which is missing from the other, symbolizing ease of dismantling.

Rascal P320 Compact

Rascal P320 Compact Specification:

Type Class	B / 2
Overall length	82.5cm / 32.5"
Overall width	58.5cm / 23"
Floor to seat height	50-52.5cm (min-max)
Seat width	43.1cm / 17"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	3.5cm / 1.4"
Maximum kerb height	4cm / 1.5"
Range*	up to 20km / 12.5 miles
The maximum safe slope	6°
Turn-around width	97cm / 3'2"
Colour	Blue

*Range on full charge and flat ground.

Main Features:

- Car transportable
- Flip-up fully adjustable armrests
- Swing-away joystick
- Transportable battery box
- Convenient joystick charging port
- Solid tyres all round
- Comfortable captains seat
- Flip-up footplate

Optional Extras:

- Optional accessory swing-away footrests
- Optional accessory elevating legrests



Rascal P320 Compact Profile:

The Rascal P320 Compact measures up to half the size of most powerchairs in today's market. With a quick release battery pack and easy to dismantle components, transportation is simply taken care of. This slim line Rear-Wheel Drive powerchair offers superb traction and outstanding manoeuvrability with a desirable turning circle of just 97cm.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E - v 3 Final-August 2014