

Rascal P312 Turnabout



*Photograph shows P312 Turnabout standard model (without seat lift).



Optional Seat Lift



Up to 29km

18 Mile
Range*



Adjustable arms



Max. Carrying
Capacity
158.8kgs / 25 Stone



Rascal P312 Turnabout

Rascal P312 Turnabout Specification:

Type Class	B / 2
Overall length	96cm / 38"
Overall width	61cm / 24"
Floor to seat height w/o seat lift	51-58.5cm (min-max)
Floor to seat height with seat lift	57-69.5cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 29km / 18 miles
The maximum safe slope	10°
Turn-around width	1.18m / 3'10"
Colour	Red

*Range on full charge and flat ground.

Main Features:

- Front-Wheel or Rear-Wheel Drive
- Flip-up width adjustable armrests
- Convenient joystick charging port
- Swing-away joystick
- High back adjustable seat
- Fold-away footplates front and rear

Optional Extras:

- Optional upgrade to a 15cm / 6" seat lift model
- Optional accessory swing-away footrests
- Optional accessory elevating legrests



Rascal P312 Turnabout Profile:

The Rascal P312 Turnabout gives flexibility to the user including a dual directional drive with a rotating seat giving you the option of Front-Wheel Drive for superb manoeuvrability indoors or Rear-Wheel Drive for extra traction outdoors. The P312 offers well defined driving positions with a swing-away joystick and a fully adjustable seat and headrest for day-long comfort.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 3 Final-August 2014