

Also Available



Matching Sofa

New Ambassador® Waterfall Back

Also Available



Static Armchair



Made in the UK



New Ambassador® Waterfall Back


The **New Ambassador®** chair is different from any other riser reclining chair with its superbly supportive seating positioning system and deep waterfall cushions providing the ultimate in seating comfort. Its unique design enables almost any individual to achieve maximum comfort, whether seated, rising or reclining.


New Ambassador® Waterfall Back Size Guide:

* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kgs	25 stone / 158.8kgs	25 stone / 158.8kgs

Main Features:

- Seminar Dual-motor Tilt-in-Space Mechanism
- Quality Craftsmanship, **BRITISH MADE** 
- Choice of sizes small, medium and large
- Available in Kilburn & Ultra Leather fabrics
- Waterfall Back Design
- New hybrid pocket sprung seat base
- New pleated wings for a timeless elegance
- Arm covers and Antimacassar as standard
- Rear lockable castors as standard
- 5 Year Guarantee



Tested by FIRA

Local Cosi Chair Retailer:

Please note images shown on this Specification Sheet are for illustrative purposes only.

The colours may vary slightly from the images shown. Riser Recliners available from STOCK in all the above fabrics.



Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset TA19 9DL

 **0800 252614**  **sales@electricmobility.co.uk**  **electricmobility.co.uk**

Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice.



MADE TO ORDER

