

Also
Available

Matching
Sofa

Ambassador

Waterfall Back



Ambassador

Waterfall Back

The Ambassador offers outstanding levels of comfort and styling with its superbly flexible seating positioning system and deep generous waterfall cushions providing the upmost essential riser recliner chair. Its unique design reaches ultimate positioning by achieving the Trendelenburg position.

Ambassador Waterfall Back Size Guide:

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	48.2cm / 19"	53.3cm / 21"	58.4cm / 23"
Seat depth	43.2cm / 17"	50.8cm / 20"	53.3cm / 21"
Arm height	16.5cm / 6½"	16.5cm / 6½"	15.2cm / 6"
Back height	66cm / 26"	68.6cm / 27"	71.1cm / 28"
Overall width	77.5cm / 30½"	81.9cm / 32¼"	83.4cm / 34"
Overall height	110.5cm / 43½"	114.3cm / 45"	119.3cm / 47"
Weight limit	21 stone / 136kgs	26.8 stone / 170kgs	26.8 stone / 170kgs

Main Features:

- Choice of Button or Waterfall Back design
- Made from Kiln Dried Canadian Maple frames
- 5 Year Guarantee on frames, mechanisms and all electrics
- Heat and massage available on request
- Multi-functional Dual-motor 6-button handset, allows independent controls of the backrest and footrest
- 2-Seater static sofa available on Made to order
- Arm covers and Antimacassar as standard
- Trendelenburg position which gives high leg elevation and raises the feet higher than the heart

Local Cosi Chair Retailer:



Leaf
Cocoa



Leaf
Oyster



Patchwork
Cocoa




Patchwork
Oyster



Seasons
Cocoa




Seasons
Rose




Seasons
Evergreen




Seasons
Oyster



Brisa
Buckskin



Brisa
Hollyhock



Brisa
Coffee Bean



Also available in
Button Back



Footrest
Extension



Brakes castors
(rear)



2-Seater
static sofa

Please note colours may vary slightly from the images shown.
Riser Recliners available from STOCK in all the above fabrics.



Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset TA19 9DL

0800 252614 sales@electricmobility.co.uk electricmobility.co.uk

