lconic



Tilt-in-Space Quad-Motor Riser Recliner





Superior & Sophisticated

- Fully adjustable powered head & lumbar support
- True footrest extension
- Longest hip to heel seat
- Fitted USB & C charging ports

Featuring Zero Gravity seating position



CONIC Quad-Motor Riser Recliner







Lasting relaxation, feel the difference

The Iconic is a Premium Riser Recliner designed by CosiChair which has a unique Tilt-in-Space (TIS) Quad-Motor function featuring the Zero Gravity position which evenly distributes the force & stress across the body, giving supported comfort from head to toe throughout movement offering greater relaxation.

Enjoy the adjustable powered head & lumbar positioning which offers personalised support for the head, neck & lower back. Benefit & take advantage of these features to watch TV, read a book or sit up to have a drink.

Designed using the finest fabrics & handmade to the highest standards, the Iconic offers first class flexibility with regards to movement & positioning, with Variable Angle Lift (VAL) adjusting the tilt angle on lift & the distinctive, automatic adjusting backrest opens up the angle for exit. This precision & smart technology allows a safe, controlled standing position with the user-friendly 9 button handset.














Features & Benefits

- Zero Gravity Seating Position
- Variable Angle Lift (VAL) to adjust the tilt angle on the lift function between 30°, 20° or 10° which reduces the seated angle when lifting & exiting the chair
- Very high leg lift giving full support behind the knee helping to relax the legs
- A 'Reset' button which has software designed to bring you back to a seated position in a safe & comfortable controlled manner from a recline or lift position
- Unique automatic adjusting backrest which helps to reduce the back angle to allow a more controlled safe standing position on exiting the chair
- Fully adjustable backrest
 Fully adjustable rise & recline headrest
 Home 'Reset' lumbar
- Longest Hip to Heel seating with an additional 3" true footrest extension
- Contemporary Lateral or Waterfall Back Design
- Industry leading 7 Year Guarantee T&C's apply*

Available to complement your home in six modern, stylish & soft fabrics. Choose from the Lateral or Waterfall sumptuous back design or indulge yourself & have both...

Other designs, upholstery & colourways to suit all tastes are available on request.







Zero Gravity

The Zero gravity position puts your body into "zero stress zone" which offers even weight distribution across the body. Your muscles relax aiding breathing & as your legs are elevated this reduces stress on the heart. The position spreads your body weight out, reducing the stress & tension on the spine, hips, knees & lower back, benefitting your health & wellbeing.



Truly improve your well-being...





^{*7} Year Guarantee mechanism & electronics warranty. Handset & transformer fair wear & tear & a lifetime wooden frame warranty. We reserve the right to amend the duration of Component Guarantee Periods identified.

Variable Angle Lift





Are you sitting comfortably?

The Iconic Riser Recliner has a Variable Angle Lift (VAL). There are 3 angles at which the chair can be set to rise at. If it isn't comfortable rising at the standard 30 degree angle, it can be set to a medium rise of 20 degrees or a "flat rise" which is just 10 degrees. The varying degrees of lift may help people that have postural problems, allowing supported & smoother transition from sitting to standing.

There is the option to configure the VAL at time of purchase.

Ask your Cosi Chair Specialist if you have specific requirements.

The Revolutionary Iconic Quad-Motor comes with the following accessories as standard:

9v battery back up transformer

Side arm mounted USB & C ports which allows you to plug in & charge your electronic devices — Twin pockets

- Antimacassar head cover

Full length arm covers

Twin handset holder for left or right-handed control
Braked castors

Choose our furniture with confidence -All riser recliner chairs are tested by FIRA - The leading authority in furniture safety by design.





Iconic Size Guide:

| | Small | Medium | Large |
|-------------------------------|-------------------|-------------------|--------------------|
| Seat height | 47cm / 18½'' | 50.8cm / 20'' | 53.3cm / 21" |
| Seat width | 48.3cm / 19'' | 50.8cm / 20'' | 53.3cm / 21′′ |
| Seat depth | 45.7cm / 18" | 50.8cm / 20'' | 55.9cm / 22'' |
| Arm height | 17.8cm / 7'' | 19.1 cm / 7½'' | 19.1 cm / 7½′′ |
| Back height | 68.6cm / 27'' | 71.1cm / 28'' | 73.7cm / 29'' |
| Hip to Heel – raised footrest | 101.6cm / 40'' | 108cm/ 42½" | 114.3cm / 45" |
| Zero Gravity position | 175.2cm / 69'' | 182.9cm/ 72" | 190.5cm / 75'' |
| Footrest extension | 7.6cm / 3" | 7.6cm / 3" | 7.6cm / 3" |
| Overall width | 76.2cm/ 30" | 81.3cm/32" | 86.4cm/ 34" |
| Overall height | 91.4cm/36" | 94cm/ 37" | 96.5cm/38" |
| Weight limit | 25 stone/ 158.8kg | 25 stone/ 158.8kg | 25 stone / 158.8kg |

Local Cosi Chair Specialist







