

The Jubilee



Tilt-in-Space Dual Motor















The Jubilee

Lateral Back

The Jubilee is our UK Premium Tilt-in-Space (TIS) Dual-motor Cosi Chair designed to maximise support throughout motion with a range of stylish back designs & upholstery options to provide exceptional postural seating positions to suit all individuals. The Jubilee is available in Waterfall or Lateral Back from stock with other back options & fabrics available Made to order.













- 2000	1 (1 to 1		(678)
Size Guide	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg

Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

Main Features & Benefits:

- Premium Seminar TIS mechanism with Variable Angle Lift* (VAL* - M & L only. Ask for more details)
- Arm covers and Antimacassar included
- 5-button power on handset
- Choice of configurable back cushions
- Heat & massage available Made to order
- Drop in seat cushions for pressure care issues available Made to order. (Ask for more details.)
- 2 or 3-Seater static sofa available on Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, BRITISH MADE







Reef Pro

Ink & Stain Resistant Technology

The unique Takumi™ technology makes Reef Pro premium fabric soft, supple, stylish & durable whilst the Mircrofoam layer guarantees the ultimate in thermal comfort.



Kilburn Mink



Kilburn **Oatmeal**



Kilburn Sky



Kilburn Plum



Kilburn Cocoa



Kilburn Zinc



Reef Pro Lionfish



Reef Pro



Reef Pro Captain Stingray

Exclusively available from your local Cosi Chair Specialist