

ELECTRIC  
MOBILITY®



# The Fenlake

 **BRITISH  
MADE**  
Quality Craftsmanship

Independent Quad-Motor



Tailored Lateral Back



Featuring Split Ottoman Extension

  
CosiChair



# The Fenlake

## Tailored Lateral Back

The Fenlake Independent Quad-Motor Riser Recliner is designed to offer exceptional comfort from head to toe. Featuring independent control, from a seated position as the legs raise, the split ottoman footrest opens up & extends to give greater & longer hip to heel support. The padded Tailored Lateral Back design with powered lumbar support & headrest, maximises the comfort level for the head, neck & lower back. The nine button handset has a 'Home Button' function which returns the chair to the standard seated position. The handset includes a practical USB port to charge a mobile phone or tablet. This Riser Recliner boasts both twin magazine & useful left & right hand control pockets.



Size Guide	Small	Medium	Large
Seat height	47cm / 18½"	50.8cm / 20"	53.3cm / 21"
Seat width	47cm / 18½"	50.8cm / 20"	54.6cm / 21½"
Seat depth	48.3cm / 19"	50.8cm / 20"	53.3cm / 21"
Arm height	16.5cm / 6½"	17.8cm / 7"	19.1cm / 7½"
Back height	71.1cm / 28"	73.7cm / 29"	76.2cm / 30"
Hip to Heel - raised footrest	94cm / 37"	97.8cm / 38½"	101.6cm / 40"
Overall width	81.3cm / 32"	83cm / 33"	88.9cm / 35"
Overall height	104cm / 41"	108cm / 42½"	111.8cm / 44"
Weight limit	25 stone / 158.8kg	25 stone / 158.8kg	25 stone / 158.8kg

### Main Features & Benefits:

- Featuring Split Ottoman Extension
- Independently powered head & lumbar support
- Upto 25 stone weight limit
- 'Reset' Home Button
- 9 Button soft touch handset complete with integral USB port
- Antimacassar & arm covers included
- Waterfall Back option available Made to Order
- 5 Year Guarantee on frames, mechanisms & all electrics
- Quality Craftsmanship, BRITISH MADE 



**Ascot Biscuit**



**Ascot Cobalt**



**Ascot Espresso**



**Ascot Mercury**



**Ascot Shark**



**Ascot Slate**