


Ultralite 480



Voted “**Best Buy**” in the Mobility Scooter Reviews 2014 by **Which?** magazine



Ultralite 480

Ultralite 480 Specification:

Type Class	A / 2
Overall length	102cm / 40"
Overall width	51cm / 20"
Floor to seat height	45-55cm (min-max)
Seat width	43.1cm / 17"
Seat depth	33cm / 13"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	4cm / 1.5"
Range*	up to 13km / 8 miles
The maximum safe slope	6°
Turn-around width	1.62m / 5'4"
Colour	Black / Red / Teal

*Range on full charge and flat ground.

Ultralite 480 Profile:

The Ultralite 480 is the ultimate transportable scooter. With its compact and lightweight design, the 480 ensures you stay independent at all times. This scooter has a proven reputation for outstanding design and reliability. Easy to transport, the 480 can be dismantled in seconds; no cables to worry about, just quick-release locks, convenient lifting handles and lightweight manageable components. Easy to stow, the 480 can be stored vertically in your home if space is an issue, and fits easily into most car boots when taken apart.



Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Charger storage in battery box
- Swivel seat
- Single heaviest component (front chassis) 13kgs / 29lbs

Optional Extras:

- A cane and crutch holder bag



Voted **“Best Buy”** by **Which?** magazine

Local Rascal Mobility Retailer: