

## Technical Note

### TRANSPORTATION OF PRODUCTS WITH SEALED LEAD ACID BATTERIES

The Rascal products below as supplied by Electric Mobility, are fitted with batteries that are Valve Regulated Lead-Acid (VRLA) batteries in Absorbent Glass Mat (AGM) design being 'Non-Spillable'

The VRLA AGM batteries supplied by Electric Mobility have passed vibration, pressure differential and free flowing acid tests under CFR 49 173.159(d) and meet **IATA Special Provisions A48 and A67.**

By the classification the batteries are declared as 'non-dangerous' and exempt from Hazardous Materials regulations and cleared for Air Travel within Passenger Aircraft.

Rascal Vie  
Rascal Vippi  
Rascal Veo Sport SR  
Rascal Vista  
Rascal Vista DX  
Rascal Pioneer  
Rascal Vecta Sport  
Rascal Vortex  
Rascal Vision  
Rascal Rio  
Rascal Razoo  
Rascal Rhythm SL  
Rascal Ryley  
Rascal Ryley SL  
Rascal Rivco  
Rascal Rivco SL  
Rascal Rueba  
Rascal Rueba CT