

# Liteway Balance Plus

Transportable scooters



Your new route to independence

**RASCAL** 

# Liteway Balance Plus

## Liteway Balance Plus Specification:

Type Class	B / 2
Overall length	107cm / 42.1"
Overall width	57.5cm / 22.6"
Floor to seat height	59cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 14Ah x 2
Upgraded battery option	12V / 22Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* standard battery pack (14Ah)	up to 14.4km / 9 miles
Range* upgraded battery pack (22Ah)	up to 26km / 16 miles
The maximum safe slope	6°
Turn-around width	1.4m / 4'5"
Colour	Blue / Graphite
Single Heaviest component	12.8kgs / 28.2lbs

\*Range on full charge and flat ground

## Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Handy storage compartment
- Tiller adjustment and can be folded down for ease of transportation
- Solid tyres all-round
- Anti-tip balance wheels
- Swivel seat



### OPTIONAL

## Optional Extras:

- Front basket
- Optional upgrade 22Ah battery pack (Range of upto 26km / 16 miles)

## Rascal Liteway Balance Plus Profile:

If you are looking for increased legroom, ease of steering and optimum manoeuvrability, the enhanced design features of the Rascal Liteway Balance Plus make this the natural choice. This lightweight, compact scooter features a two-piece chassis designed to dismantle with ease and a removable transaxle, all of which fit comfortably into the boot of a car. Available with an optional battery upgrade to 22Ah, that can give up to a 16 mile range.

Local Rascal Mobility Retailer:

