

# Liteway 4 Plus

Transportable scooters



Your new route to independence

**RASCAL** 

# Liteway 4 Plus

## Liteway 4 Plus Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	57.5cm / 22.6"
Floor to seat height	59cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* standard battery pack (14Ah)	up to 14.4km / 9 miles
Range* upgraded battery pack (22Ah)	up to 26km / 16 miles
The maximum safe slope	6°
Turn-around width	1.6m / 5'2"
Colour	Blue / Graphite
Single Heaviest component	14.9kgs / 32.8lbs

\*Range on full charge and flat ground

## Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Handy storage compartment
- Tiller adjustment and can be folded down for ease of transportation
- Solid tyres all-round
- Swivel seat

### OPTIONAL



### Optional:

- Optional accessory a front basket
- Optional upgrade 22Ah battery pack

## Rascal Liteway 4 Plus Profile:

This high performance all-rounder is daringly different thanks to the unique styling of the lightweight tubular frame and innovative latch-lock system on the two-piece chassis. Thoughtful engineering makes it so easy to dismantle and stow into an average sized boot, leaving plenty of room for other necessities - smaller, individual sections are easier to lift and provide even greater flexibility for when storing away purposes are essential. Available with an optional battery upgrade to 22Ah, that can give up to a 16 mile range.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E - v 5 April 2019