

Liteway 8 MaxGrip

Transportable scooters



Your new route to independence

RASCAL 

Liteway 8 MaxGrip

Liteway 8 Specification:

Type Class	C / 3
Overall length	126cm / 49.6"
Overall width	60cm / 23.6"
Floor to seat height	63cm (min)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	3 sets of 12V 14Ah x 2
Maximum speed	12.8km/h / 8mph
Ground clearance	8.5cm / 3.3"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 38km / 23.8 miles
The maximum safe slope	10°
Turn-around width	1.98m / 6'5"
Colour	Blue / Graphite
Single Heaviest component	18.1kgs / 40lbs

*Range on full charge and flat ground

Main Features:

- Rear coil spring suspension
- High efficiency LED lighting
- Lightweight aluminium frame
- Fold-down tiller for easy storage & transportation
- Luxury seating with swivel function & flip up / width adjustable armrests

OPTIONAL ACCESSORY



Optional Accessories:

- Front basket
- Super compact seat
- Off-board charging station

Rascal Liteway 8 Max Grip Profile:

At last, a Class 3, transportable, aluminium scooter that has been uniquely designed to go the extra mile. The Rascal Liteway 8 Max Grip provides the versatility of a transportable scooter with the performance & range of a full Class 3 (road legal) scooter normally too heavy & cumbersome to transport.



Local Rascal Mobility Retailer:

• **Max Grip LSD - Limited Slip Differential Gearbox Features:**

- Improved grip, handling & safety
- Unique take-apart construction & fully road legal
- Improved energy efficiency - increased range