

# 388 S

Mid size scooters



Your new route to independence

**RASCAL** 

# 388 S

## 388 S Specification:

Type Class	B / 2
Overall length	125cm / 49"
Overall width	54cm / 21"
Floor to seat height	59-64cm (min to max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 33Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	10cm / 4"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	1.7m / 5'6"
Colour	Red / Silver

\*Range on full charge and flat ground

## Main Features:

- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable height adjustable, sliding swivel seat
- Narrow compact chassis

### UPGRADE



## Upgrade:

- The familiar yet distinct Rascal 388 XL, is part of this family of reliable scooters. With a maximum speed of 6 mph and a ground clearance of 13cm, enquire today about the Class 3 option.

## Rascal 388 S Profile:

The Rascal 388 S (Standard) is a well established, great value for money mid-size scooter to suit your everyday mobility needs. At 21" wide, it is the narrowest width of any scooter in its class & is equally at home inside a busy supermarket as it is on the pavement. With a comfortable swivel seat, smooth controls & adjustable tiller, this high performance, reliable scooter is your ideal mobility partner.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E- v7 August 2021

 [electricmobility.co.uk](http://electricmobility.co.uk)