

# 388 XL

Mid size scooters



Your new route to independence

**RASCAL**



# 388 XL

## 388 XL Specification:

Type Class	C / 3
Overall length	125cm / 49"
Overall width	54cm / 21"
Floor to seat height	61 - 65cm (min-max)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage & capacity	12V / 33Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	13cm / 5"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	1.7m / 5'6"
Colour	Black / Red

\*Range on full charge and flat ground

## Rascal 388 XL Profile:

The Rascal 388 XL is the perfect Class 3 scooter with 6 mph top speed, stylish 11" wheels, full road lighting, a wing mirror and with its narrow chassis it combines safety with performance and manoeuvrability. Sturdier in build and ideal for people requiring additional comfort on those slightly longer journeys - a superb scooter for days out and leisure trips. With a comfortable high back seat with headrest, width adjustable armrests and tiller, this scooter can be 'adapted' to suit the individual's requirements.

Local Rascal Mobility Retailer:

## Main Features:

- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis
- Lift-up tiller head for easy access
- Captain's seat
- Smooth controls
- One wing mirror as standard
- Front basket as standard

DISCOVER



## Other Models:

- Enquire about the Rascal 388 S today - A narrow Class 2, 4mph scooter for the everyday.



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E. - v 7 August 2021