

Rascal 850

Performance scooters



Your new route to independence

RASCAL 

Rascal 850

850 Specification:

Type Class	C / 3
Overall length	127cm / 50"
Overall width	65cm / 26"
Floor to seat height	62cm
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	15cm / 6"
Maximum kerb height	10cm / 4"
Range*	up to 38km / 24 miles
The maximum safe slope	9°
Turn-around width	1.70m / 5'6"
Colour	Black / Red

*Range on full charge and flat ground

Rascal 850 Profile:

Whether driving at 4 mph on pavements or 8 mph on roads, the Rascal 850 brings together compact design with Class 3 road legal specification. Available in black / red and with 50Ah batteries as standard you can travel in confidence and style. A fully adjustable sculpted sliding seat with added lumbar support combined with full suspension ensures a smooth ride, whatever the terrain.

Local Rascal Mobility Retailer:

Main Features:

- Five-position flip-up tiller
- Flip-up width adjustable armrests
- All round coil suspension
- Height adjustable headrest
- Backrest angle alteration
- Sturdy front bumper
- One wing mirror as standard

OPTIONAL ACCESSORY

Optional Extras:

- Rascal armrest bag
- Cane & Crutch holder bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v 8b May 2021