

P320 Compact

Compact powerchairs



Your new route to independence

RASCAL 

P320 Compact

P320 Compact Specification:

Type Class	B / 2
Overall length	82.5cm / 32.5"
Overall width	58.5cm / 23"
Floor to seat height	50-52.5cm (min - max)
Seat width	43.1cm / 17"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	3.5cm / 1.4"
Maximum kerb height	4cm / 1.5"
Range* - standard battery pack	up to 20km / 12.5 miles
The maximum safe slope	6°
Turn-around width	97cm / 3'2"
Colour	Blue
Single Heaviest component	27.5kg / 60.6lbs

*Range on full charge and flat ground

Main Features:

- Car transportable
- Flip-up fully adjustable armrests
- Swing-away joystick
- Transportable battery box
- Convenient joystick charging port
- Solid tyres all round
- Comfortable Captain's seat
- Flip-up footplate



OPTIONAL ACCESSORY

Optional Extras:

- Armrest bag
- Oxygen bottle bag
- Cane & Crutch holder bag
- Swing-away or elevating legrests also available

Rascal P320 Compact Profile:

The Rascal P320 Compact measures up to half the size of most powerchairs in today's market. With a quick release battery pack and easy to dismantle components, transportation is simply taken care of.

This slim line Rear-Wheel Drive powerchair offers superb traction and outstanding manoeuvrability with a desirable turning circle of just 97cm.



Local Rascal Mobility Retailer:



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E - v5 Final-March 2019