

Your new route to independence

P327 Mini

P327 Mini Specification:

Type Class	B/2
Overall length	109cm / 42.9"
Overall width	62.4cm/24.5"
Floor to seat height	52cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs/21 stone
Battery voltage and capacity	12V/34Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm/2"
Maximum kerb height	5cm/2"
Range* - standard battery pack	up to 35km / 22 miles
The maximum safe slope	6°
Turn-around width	1.18m/3'8"
Colour	Red / Teal
Single Heaviest component	62kg/136.7lbs

^{*}Range on full charge and flat ground

Main Features:

- Outstanding manoeuvrability
- Flip-up fully adjustable armrests
- Full suspension
- Convenient joystick charging port
- Dual density, puncture resistant drive wheels for maximum comfort on the
- Swing-away joystick
- High back adjustable Captain's seat



Optional Extras:

- Armrest bag
- Oxygen bottle bag
- Cane & Crutch holder bag



Rascal P327 Mini Profile:

The Rascal P327 Mini is a compact 4 mph Mid-Wheel Drive indoor/outdoor powerchair, with full suspension and comfortable scooter style seating. The P327 Mini is the ideal choice for indoor and outdoor use where space is restricted. It delivers great performance with its 34Ah batteries contributing to its outstanding manoeuvrability and a smooth ride.

Discover the P327 Mini Seat Lift option with a column lift raising you by 7.5 inches which delivers a smooth lift action.

Local Rascal Mobility Retailer:







