

P327

Mid-wheel powerchairs



Your new route to independence

RASCAL 

P327

P327 Specification:

Type Class	B / 2
Overall length	112cm / 44"
Overall width	69cm / 27"
Floor to seat height	54- 56.5cm (min-max)
Seat width	47cm / 18.5"
Seat depth	40.6cm / 16"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20m / 4'
Colour	Red / Teal
Single Heaviest component	

*Range on full charge and flat ground

Rascal P327 Profile:

The Rascal P327 is a 4 mph Mid-Wheel Drive outdoor/indoor powerchair, with durable full suspension and comfortable cushioned Captain's seat. The P327 delivers great performance with its 55Ah batteries contributing to its excellent range, outstanding manoeuvrability and a smooth ride.

Enquire today about upgrading to a larger seat at 20.5". P327 Mini, P327 XL and Seat Lift models also available.

Local Rascal Mobility Retailer:


Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable captains seat
- Flip-up footplate

OPTIONAL ACCESSORY

Optional Extras:

- Armrest bag
- Oxygen bottle bag
- Cane & Crutch holder bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E. - v 5 Final-May 2019