

# Rascal P200

Mid-range powerchairs



Your new route to independence

**RASCAL** 

# Rascal P200

## P200 Specification:

Type Class	B / 2
Overall length	109cm / 42.9"
Overall width	65.1cm / 25.6" - 73.7cm / 29" (min - max)
Floor to seat height	57cm (min)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 26km / 16.2 miles
The maximum safe slope	8°
Turn-around width	1.32m / 4'3"
Colour	Silver

\*Range on full charge and flat ground

## Rascal P200 Profile:

The Rascal P200 4 mph Rear-Wheel Drive powerchair suitable for indoor and outdoor use comes with 55Ah batteries as standard and features a flexible seating system with a removable seat cushion enabling bespoke adaptations. The capacity to adjust the seat tilt, backrest angle and the horizontal strap tension for total comfort achieves the ideal seating position.

Local Rascal Mobility Retailer:

## Main Features:

- Adjustable upholstery on back
- Removable armrests
- Coil spring suspension
- Convenient joystick charging port
- Solid tyres all round
- Manual angle adjustable seat tilt
- Swing-away legrests

### OPTIONAL ACCESSORY



## Optional Extras:

- Optional accessory elevating legrests
- Option for Kerb Climbers
- Option of an Attendant control bracket



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E - v 5 Final-April 2019