

# Veo Sport

Transportable scooters



Your new route to independence

**RASCAL** 

# Vevo Sport

## Vevo Sport Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	47-49cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	129.2kgs / 20 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.3"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 21km / 13.5 miles
The maximum safe slope	8°
Turn-around width	1.40m / 4'6"
Colour	Aquamarine / Platinum Grey
Single Heaviest component	20.2kgs / 44.5lbs

\*Range on full charge and flat ground

## Rascal Vevo Sport Profile:

It's all about the drive. The Rascal Vevo Sport is fully loaded with features. A modern, sporty design with class leading all-round suspension. With 22Ah batteries, a superb capacity for a scooter of its size, affording a range of up to 13.5 miles. A padded delta tiller and luxury seat, detailed with unique white accent stitching, the Vevo Sport is a real pleasure to ride.

The Vevo Sport's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension, this is the travel scooter with the edge.

Local Rascal Mobility Retailer:

## Main Features:

- Car transportable
- All-round suspension
- Stylish sports seat
- Dual compound puncture proof tyres
- Comfortable delta tiller
- Rear safety reflectors
- Colours in latest high-tech, two-tone matt finishes

OPTIONAL



## Optional:

- Front lighting system
- Now available with lightweight 14Ah battery pack (Range of up to 18.5 km / 11.5 miles)
- Pneumatic tyres



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E - v 4 February 2019