



The **DIETZ**[®] TAIiMA Rollator Range

What a rollator should be



TAIMA M

Optimal stability,
foldability and
modern looks

Key Features:

- Walking stick holder
- String shopping bag (up to 5kg)
- Tilting aid to overcome obstacles
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience



TAIMA XC

Hardwearing, rugged design

Key Features:

- Walking stick holder
- String shopping bag (up to 5kg)
- Tilting aid to overcome obstacles
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience



TAIMA XC images shown with optional backrest (not included as standard)

TAiMA M-GT

Fully featured flagship model

Key Features:

- Lightweight - 5.5kg* / 5.8kg**
- Back rest
- Walking stick holder
- String shopping bag (up to 5kg)
- Tilting aid to overcome obstacles
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience
- Soft polyurethane tyres for smoother wheel movement



TAiMA S-GT

Smaller in size, fully featured flagship model

Key Features:

- Lightweight - 5.3kg* / 5.6kg**
- Back rest
- Walking stick holder
- String shopping bag (up to 5kg)
- Tilting aid to overcome obstacles
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience
- Soft polyurethane tyres for smoother wheel movement



Specifications



Product Name	TaiMA M	TaiMA GT		TaiMA XC
Size	Standard	S-GT	M-GT	Standard
Colour	Metallic Brown	Metallic Sepia	Metallic Reseda Green	Metallic Red
Overall height (cm)	79 - 95	76 - 86	79 - 95	85 - 95
Height folded (cm)	81	76	81	88
Seat height (cm)	59	51	59	61
Seat width (cm)	45	39	45	45
Seat depth (cm)	20	20	20	23
Overall length (cm)	68	62	68	75
Overall width (cm)	63	58	63	67
Width folded (cm)	23	23	23	27
Max. Capacity (kg)	150	150	150	150
Push handles adjustable	5 steps / 4 cm	5 steps / 2.5 cm	5 Steps / 4 cm	5 Steps / 2.5 cm
Overall weight (kg)	6.3	5.3* / 5.6**	5.5* / 5.8**	9.1**

* Excludes string shopping bag, back rest and walking stick holder. ** Including string shopping bag, back rest and walking stick holder.

Accessories



Rollator Storage Bag

Made of durable, black, water-proof nylon, washable up to 30°C, easy to open and close with one hand. The bag can be folded away to save space.



Comfort Backrest

Upholstered, padded backrest, for added comfort. Washable and easily fitted in place of the standard backrest.



Quadro X Handlebar Light

Super bright, battery powered LED light. Clips easily to either handlebar, single button operation - see and be seen in the dark.

To place an order, please contact our Customer Service team:

Call: 0800 252 614

email: sales@electricmobility.co.uk

**ELECTRIC
MOBILITY**



Electric Mobility Euro Limited, Canal Way, Ilminster, Somerset TA19 9DL
Registered in England Number 2419231

Issue 1 May 2017 | All rights reserved E.&O.E