

P327 Mini Seat Lift

Mid-wheel powerchairs



Your new route to independence

RASCAL 

P327 Mini Seat Lift

P327 Mini SL Specification:

Type Class	B / 2
Overall length	110cm / 43.6"
Overall width	62.4cm / 24.5"
Floor to seat height Lowest to highest elevation position	56-75cm (min-max)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 34Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 35km / 22 miles
The maximum safe slope	6°
Turn-around width	1.19m / 3'9"
Colour	Red / Teal
Single Heaviest component	65kg / 143.3lbs

*Range on full charge and flat ground

Main Features:

- 7.5 inch column lift
- Full suspension
- Outstanding manoeuvrability
- Flip-up fully adjustable armrests
- Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Dual density, puncture resistant drive wheels for maximum comfort on the move
- Maximum floor to seat height when fully elevated - 75 cm / 29.5"



Rascal P327 Mini Seat Lift Profile:

The Rascal P327 Mini Seat Lift is a medium-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, where space is restricted & with durable full suspension. The P327 Mini Seat Lift has a range of up to 22 miles and delivers great performance with 34Ah batteries. This model includes a comfortable Captain's seat with a 7.5" column lift which delivers a smooth lift action and reaches a maximum floor to seat height of 75 cm when fully elevated.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E - v 4 Final-May 2019