

# Rascal P327 XL Seat Lift



Max. Carrying  
Capacity  
182kg / 28 Stone

**RASCAL**  
MOBILITY

**Mid-wheel powerchairs**

**Your new route to independence**

# Rascal P327 XL Seat Lift

## Rascal P327 XL Seat Lift Specification:

Type Class	B / 2
Overall length	112cm / 44"
Overall width	69cm / 27"
Floor to seat height	62-87cm (min-max)
Lowest to highest elevation position	
Seat width	52cm / 20.5"
Seat depth	45.7cm / 18"
Maximum carrying capacity	182kgs / 28 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20m / 4'
Colour	Red / Teal

\*Range on full charge and flat ground.

## Main Features:

- 10 inch scissor lift mechanism
- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Maximum floor to seat height when fully elevated - 87 cm / 34"

## Optional Upgrade

- Larger seat available - 56cm / 22"



## Rascal P327 XL Seat Lift Profile:

The Rascal P327 XL Seat Lift is a heavy-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, with durable full suspension & outstanding manoeuvrability. The P327 XL Seat Lift has a range of up to 20 miles and delivers great performance with 55Ah batteries. This model includes a comfortable Captain's seat with a 10" scissor lift mechanism. The scissor lift delivers smooth and quiet lift action and reaches a maximum floor to seat height of 87 cm when fully elevated.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E - v 4 Final - July 2018