


The

Banwell

Tilt-in-Space Dual-motor

Made
in the UK



Banwell

Waterfall Back

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



Riser Recliners available from STOCK in the following fabric



Colours may vary slightly
from images shown which
are for illustrative purposes only.



The Banwell is a powerful & significant Riser Recliner chair for the heavier or wider customer that provides the upmost reassurance & comfort.

40 stone Banwell with the 25 stone leg lift – option 1

This has been designed with the larger customer in mind, & those who may require a footboard suitable to support extra carrying capacity. Unlike other standard riser recliner chairs available, which can only carry approximately 20% total lifting capacity on the footboard, the CosiChair Banwell carries up to 25 stone on its footboard area. It is a Dual-motor chair controlling the base and back with a further 2 motors powering the leg lift only.

35 stone Banwell with Tilt-in-space – option 2

This Dual motor recliner is for the heavier customer who requires unparalleled comfort & support of the head and back area. As you recline the chair, the seat to back angle doesn't alter allowing for full support throughout the chairs motion, it then allows the back to be separately adjusted for the ultimate comfort and positioning.

*Available to order in Seasons Cocoa. Other fabric options available on a Made to order basis - please ask for more details.

Main Features:

- 2 variations of mechanism available:
35 stone Tilt-in-Space or 40 stone with 25 stone leg lift
- Arm covers and Antimacassar included
- Waterfall Back design with wings
- High Density Foam seat base
- Pressure care cushions available on request*
- Made to order & Made to measure available
- 2 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE**

Banwell Size Guide for stocked product:

	40 stone option 1	35 stone option 2
Seat height	48.2cm / 19"	48.2cm / 19"
Seat width	64.8cm / 25½"	64.8cm / 25½"
Seat depth	53.3cm / 21"	53.3cm / 21"
Arm height	20.3cm / 8"	20.3cm / 8"
Back height	88.9cm / 35"	88.9cm / 35"
Overall width	81.3cm / 32"	81.3cm / 32"
Overall height	122cm / 48"	122cm / 48"
Weight limit	40 st / 254kg	35 st / 222kg

Local Cosi Chair Retailer:

