

The
Jubilee

Tilt-in-Space Dual-motor

Made
in the UK






Jubilee

Lateral Back

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



Riser Recliners available from STOCK in the following fabrics



Colours may vary slightly from images shown which are for illustrative purposes only.

The Jubilee is our UK Premium Tilt-in-Space (TIS) Dual-motor Cosi Chair designed to maximise support throughout motion with a range of stylish back designs & upholstery options to provide exceptional postural seating positions to suit all individuals. The Jubilee is available in Waterfall or Lateral Back from stock with other back options & fabrics available Made to order.

Local Cosi Chair Retailer:

Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift* (VAL* - M & L only. Ask for more details)
- Arm covers and Antimacassar included
- 5-button power on handset
- Choice of configurable back cushions
- Heat & massage available Made to order
- Drop in seat cushions for pressure care issues available Made to order. (Ask for more details.)
- 2 or 3-Seater static sofa available Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE**

Jubilee Lateral Size Guide:

* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg

