

The

Kensey

Tilt-in-Space Single-motor

Made
in the UK



Available



Made to
Measure

**MADE
TO
ORDER**




CosiChair

Kensey

MADE TO ORDER

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by design.



The Kensey is available Made to order or Made to measure ONLY.



Lateral Back.
2 or 3-Seater static sofa available on Made to order

Colours may vary slightly from images shown which are for illustrative purposes only.

The Kensey is our UK Premium Tilt-in-Space (TIS) Single-motor Cosi Chair providing a smooth rise and recline movement. The choice of 4 sumptuous Back designs deliver support to the user's back, shoulders and neck areas throughout travel. Choose from Button, Lateral, Mini Lateral or Waterfall Back options.

Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift* (VAL* - M & L only. Ask for more details)
- Arm covers and Antimacassar included
- 3-button power on handset
- Choice of configurable back cushions
- Heat & massage available Made to order
- Drop in seat cushions for pressure care issues available Made to order. (Ask for more details.)
- Choice of any back design or style of fabric on a Made to order basis **only**
- 2 or 3-Seater static sofa available Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE**

Local Cosi Chair Retailer:

Kensey Size Guide:

* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order only

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	66cm / 26"	68.6cm / 27"	71.1cm / 28"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg



MADE TO ORDER

