

The

Walden

Wallhugger Single-motor

Made
in the UK



Walden

Waterfall Back

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



Riser Recliners available from STOCK in the following fabrics



Colours may vary slightly from images shown which are for illustrative purposes only.

Main Features:

The Walden Cosi Chair is a relaxing single motor wallhugger chair with a pocket sprung seat base and a deep fibre filled Waterfall backrest. The key feature is it can sit very close to a wall or solid object (with a distance of 4") and still offer a full recline, it manages this by the whole wooden frame moving over its own axis. This chair is ideal for small rooms or confined areas where space is limited.

- A 3 button power on hand control
- Waterfall Back design
- Pocket sprung seat base
- Scroll arms with wings
- Arm covers and Antimacassar included
- Other fabrics are available Made to order
- 2 or 3-Seater static sofa available Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE**

Walden Size Guide:

	One size
Seat height	50.8cm / 20"
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Arm height	19.1cm / 7½"
Back height	68.6cm / 27"
Overall width	83.8cm / 33"
Overall height	110.5cm / 43½"
Weight limit	25 stone / 158.8kg

Local Cosi Chair Retailer: