

Smilie

Transportable folding scooter



Your new route to independence

RASCAL 

Rascal **Smilie** Profile:

The **Rascal Smilie**'s easy-to-fold mechanism allows you to fold the scooter manually for easy transportation and storage in just a few simple steps.

Featuring four equally sized 8" wheels at each corner for extra stability, the **Smilie** is the first folding scooter with a delta-bar tiller. The compact & clever styling makes it easy to manoeuvre, lightweight and perfect for travelling at home or overseas.

Enjoy the convenience & impressive durability...



**...at home
or away!**



Your perfect travel companion

Main Features:

- Car, cruise, plane or train transportable
- Lightweight & removable Lithium-ion Battery technology
- Airline Compliant 11.6Ah battery weighs just 2.8kgs / 6.2lbs
- Long Range 17.4Ah battery weighs 3.4kgs / 7.5lbs
- Fold down backrest & easy to remove adjustable armrests
- Complete aluminium design
- Delta-bar tiller
- 4 x 8" solid tyres all-round with impressive stability
- 14.5" seat height that offers real comfort & support for long journeys
- Folded Dimensions (LxWxH) 77cm x 50cm x 51cm / 30.3" x 19.7" x 20"



Clever, practical **front handle** increases accessibility

Boot Transportable



Air travel & CAA Lithium-ion Battery Legislation

Airline Compliant 11.6Ah battery weighs just 2.8kgs / 6.2lbs for overseas travel & **can be taken on a plane.**

The Long Range 17.4Ah battery cannot be taken on a plane.

Choose one of each, so you can travel up to 11.2 miles at home or 7.5 miles when you are away with the Airline Compliant battery.



Foldable & compact, simple to store or transport.
Shown here in Horizontal Fold Position & Vertical Storage Position.

Smilie

Specification:	Smilie 11.6Ah	Smilie 17.4Ah
Type Class	B / 2	✓
Overall length	96.2 cm / 37.9"	✓
Overall width	50cm / 19.7"	✓
Overall height	96cm / 37.8"	✓
Floor to seat height	56cm / 22"	✓
Seat width	41.9cm / 16"	✓
Seat depth	34cm / 13.4"	✓
Folded Dimensions LxWxH	77cm / 50cm / 51cm (30.3" / 19.7" / 20")	✓
Maximum carrying capacity	115kgs / 18 stone	✓
Battery voltage and capacity	11.6Ah x 1	17.4Ah x 1
Maximum speed	6.4km/h / 4mph	✓
Ground clearance	4.5cm / 1.8"	✓
Maximum kerb height	5cm / 2"	✓
Range*	up to 12km / 7.5 miles	up to 18km / 11.2 miles
The maximum safe slope	6°	✓
Turn-around width	1.35m / 4'4"	✓
Colour	Peach / Sky Blue	Peach / Sky Blue
Product weight (Including battery**)	29.9kgs / 65.9 lbs	30.6kgs / 67.5 lbs

*Range on full charge and flat ground

** Product weight excludes armrests

Local Rascal Retailer:



Specifications subject to change without prior notice
Images are for illustrative purposes only
SM ONLY E. and O.E- v 3 b August 2019