




Vortex Technical Specification



Authorised Electric Mobility
Retailer:

Type Class	C / 3
Overall length	146 cm / 57.5"
Overall width	65 cm / 25.6"
Type of Seat	Captain (enhanced) - swivel, head rest, reclining/folding back & arms with width adjustment and seat- slide function.
Maximum carrying capacity	180kg / 28 stone
Battery voltage and capacity	12V / 75Ah x 2
Maximum speed	12.87km/h / 8mph
Ground clearance	14 cm / 5.5"
Maximum safe climbing kerb height	10 cm / 4"
Maximum safe descend-able kerb height	10 cm / 4"
Range*	Up to 41km / 25.5 miles*
The maximum safe slope	10°
Turn-around width	2.2m / 7'2"
Turning Diameter	4.0m / 13'01"
Lighting	LED, Front and Rear
Wheel Diameter	34 cm / 13.4"
Type of tyres	Run-flat pneumatic RunOn tyres
Colour	Titanium Blue / Metallic Heather

* Range on full charge and on flat ground.

ELECTRIC
MOBILITY

Product images are for illustrative purposes only and the colour may differ slightly on the actual product due to the limitations of photography & the printing process. Specifications subject to change without prior notice.

E. and O.E- V2 Aug 2021



Vortex

Take your journeys to the next level



The Vortex is a revolutionary, unique & dynamic 8mph scooter ready to take your mobility & lifestyle to the next level. Combining new edge design with the latest technology, the Vortex is built to perform with a host of user-friendly features.

Style, substance & safety

Plenty of legroom coupled with full support & comfort from its Captain's Seat with sleek white stitching & modern design.

Performance suspension ensures a driving experience that is more dynamic & comfortable than ever before for prolonged journeys of upto 25.5 miles* over all terrains.

Super bright LED lights & extra mirrored running light for increased visibility. 2 rear view mirrors & USB charging point offers the ultimate in peace of mind motoring.

Enjoy the drive

Max Grip limited slip differential gear box & 900W high torque motor enable the Vortex to keep going where other scooters would struggle. The seat & tiller are adjustable to further enhance the ergonomics for lumbar support.

Engineered for manoeuvrability & usability with uncompromising features & stylish, versatile design.



Vortex Driving Characteristics

- Robust & sturdy
- **Max Grip** limited slip differential gearbox
- Run flat **pneumatic RunOn** tyres
- upto **28 stone** max carrying capacity
- Bright **LED** headlights
- Dedicated waterproof superfast **USB port** for charging your mobile phone
- Two striking colours
Titanium Blue or **Metallic Heather**



Everyone has the right to live as full and independent a life as possible.

Due to a policy of continual improvement, Electric Mobility Euro Ltd reserves the right to change product specifications without prior notice. Whilst every care is taken to ensure that the information included in this brochure is correct, no legal responsibility can be accepted for any inaccuracy. Images are for illustrative purposes only. Details are correct at the time of going to print. Full Terms and Conditions are available from your Authorised Electric Mobility Retailer or online at electricmobility.co.uk