

ELECTRIC MOBILITY



Authorised Electric Mobility Retailer:

Due to a policy of continual improvement, Electric Mobility Euro Ltd reserves the right to change product specifications without prior notice. Whilst every care is taken to ensure that the information included in this brochure is correct, no legal responsibility can be accepted for any inaccuracy. Images are for illustrative purposes only. Details are correct at the time of going to print. Full Terms and Conditions are available from your Authorised Electric Mobility Retailer or online at electricmobility.co.uk

Product images are for illustrative purposes only and the colour may differ slightly on the actual product due to the limitations of photography & the printing process. Specifications subject to change without prior notice. E. and O.E- V4 Final October 2025



vecta sport



VECTA

Sharp shape and rock-solid build quality make this a modern classic.

Welcome to the Vecta Sport, a revolutionary scooter that combines new edge design with the latest technology to make a scooter that looks good and performs even better. The result of three years' development, the Vecta Sport is built to perform with a host of user-friendly features, like the 25-stone carrying capacity, the "RunOn" run-flat tyres and revolutionary "Max Grip" transmission, to keep you going when other scooters cannot. Wherever you are, with a range of up to 27 miles on a single charge, there is nothing to stop you.

VECTABILITY





DESIRABILITY



Style comes as standard.

Designed without compromise, the Vecta Sport combines echoes of iconic styling with a modern clean silhouette. But it doesn't just look good, the Vecta Sport is practical too. Available in two contemporary colours, Atomic Red and Cobalt Blue, the Vecta Sport is in a class of its own.

CAPABILITY



Travel first class, every time.

The Vecta Sport offers the best of both worlds. Robust and sturdy on the road, agile enough for the pavement and at only 52 inches long, it doesn't require a large storage space. Whilst the new high power, high torque 600 Watt motor will tackle even the steepest hills with ease.



Take a seat and relax in style.

Sit behind the adjustable tiller and you will benefit from a protected position that offers an enhanced sense of safety.

Behind the tiller you will find a useful storage pocket to hold small items and a dedicated USB port to charge your mobile phone.

Wherever you go, you will find the unique mirrored running light, which is positioned high on the tiller providing style and aiding visibility. This combined with super bright LED headlights ensures you can see, and more importantly, be seen at all times.

USABILITY





MANOEUVRABILITY

Born to run and run and run.

The Vecta Sport has performance and substance. The 600 Watt high power motor gives you a range of up to 27 miles, perfect for those long days out visiting friends and family or exploring the countryside.

The powerful high torque, super efficient motor coupled with the latest "Max Grip" transmission, featuring a limited slip

differential gearbox, ensures that the Vecta Sport keeps going where most scooters would struggle.

And in the event of a flat? You just keep going. The Vecta Sport offers the ultimate in peace of mind motoring with its "RunOn" tyres ensuring you get home safely, even if you get a puncture.

Alternative Tyre option:
Puncture proof infill tyres available.
Ask for more details.



Enjoy the ride, wherever you go.

With plenty of leg room and 3-way adjustable seating, the Vecta Sport can be tailored to suit every shape and size.



ADJUSTABILITY




Add to that, a highly-padded seat with lumbar support and flip up armrests for easy access and you will always travel in total comfort.

The Vecta Sport is completely versatile.

The seat has many different positions, moving forwards and backwards, up and down, and the tiller is also adjustable, tilting towards and away from the rider to further enhance the ergonomics and your sense of safety and security.


Vecta Sport Technical Specification.



RUN-ON TYRES




SUSPENSION




600 WATT HIGH POWER MOTOR



MAX GRIP - LIMITED SLIP DIFFERENTIAL GEARBOX



USB PORT



LED



Vecta Sport meets the requirements as a Class 3 Scooter with necessary user equipment as defined with the Invalid Carriage Act (Amended 2025)

Type Class	C / 3
Overall length	132 cm / 52"
Overall width	60 cm / 23.6"
Floor to seat height	62 - 67 cm (min - max)
Type of Seat	Captain (enhanced) - swivel, head rest, reclining/folding back & arms with height & width adjustment & seat- slide function.
Maximum carrying capacity	160kg / 25 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	12.87km/h / 8mph
Ground clearance	10 cm / 4"
Maximum safe climbing kerb height	10 cm / 4"
Maximum safe descend-able kerb height	10 cm / 4"
Range*	Up to 44km / 27.3 miles
The maximum safe slope	10°
Turn-around width	1.7m / 5'5"
Turning Diameter	3.47m / 11'38"
Lighting	LED, Front and Rear
Wheel Diameter	30.5 cm / 12"
Type of tyres	pneumatic - "RunOn" Tyres
Colour	Atomic Red / Cobalt Blue

*Range on full charge and on flat ground.

Every scooter carries our reputation.

Electric Mobility is a UK company with a proven track record stretching back over 35 years. In that time, we have been developing and supplying products that help people with their mobility and daily living.

Today, all that experience has gone into the design and manufacture of the Vecta Sport, leading the way and setting a new benchmark in scooter technology.

Electric Mobility believes everyone has the right to live as full and independent a life as possible. The Vecta Sport has been designed with this in mind.

Vecta Sport Dimensions.

