

MC2RS Durable Rear-Wheel Drive Powerchair



Next level in style

EXCEPTIONALLY SPORTY DESIGN

MC2 RS: Power and good looks together at last. This special edition of the beautifully engineered MC2 takes you to the next level in style and power, with an MC2 RS exclusive seating system, 15" black anodised drive rims, optional wide tyres for outdoors as well as customisation and accessories galore!






Stylish, sporty & powerful with red & black colour scheme

Highly manoeuvrable for indoor & outdoor use

outdoor use Versatile-can be easily customised for

your wants & needs



High quality German design & manufacture

Crash tested—Complies with international ISO crash test requirements

*Range: Dependent on weight & environment. Full charge & on flat ground



15" drive wheels



Side guard lighting with custom engraving



Black- Red RS design



Recaro seat in black leatherette - optional

MC2 RS Specification

Seat width	380 – 650 mm
Seat depth	400 – 560 mm
Backrest height	530 – 640 mm
Legrest length	350 – 500 mm
Armrest height	180 – 290 mm
Seat height - front	530– 630 mm, (depending on tyre combinations)
Overall length – with legrests	Approx. 1130 mm
Overall length – without legrests	Approx. 830 mm
Overall width	635 mm
Wheel size front	260 x 70mm
Wheel size rear	380 x 75mm
Maximum obstacle climbing height	Approx. 60mm
Data subject to changes	

Turning radius	Approx. 840 mm
Driving Speeds	4mph / 8mph
Driving range	22 - 25 miles / 35 – 40 kms with 73 Ah (20h) batteries ¹
Chair weight	Approx. 100 kgs, no batteries
Maximum client weight	160 kgs / 25 stone
Motor power rating	350 watt 4 pole German motors
Maximum gradient	15°
Mechanical seat tilt	0° to 10°
Seat raise	300mm
Power seat tilt	0° to 50°, tilt in space
Mechanical backrest	-10° to + 30°
Power backrest recline	-10° to + 50°
1 Distance will be affected dimensions in mm speeds in mph (km/h) by client weight, terrain weights in kg dimension tolerances ± 10 mm and battery condition wheels in inches ± 2°	



MEYRA We reserve the right to make technical changes and cosmetic changes to our products. These changes may not be reflected in this brochure. T's & C's as stated at electricmobility.co.uk/pages/terms-amp-conditions/53/ shall apply