

Vantage X

Mid-size scooters



Your new route to independence



Vantage X

Vantage X Specification:

Type Class	C / 3
Overall length	120cm / 47.2"
Overall width	58.5cm / 23"
Floor to seat height	53-58cm
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity**	149.7kgs / 23.5 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	6.5cm / 2.6"
Maximum kerb height	7.6cm / 3"
Range*	up to 42km / 26 miles
The maximum safe slope	7°
Turn-around width	1.40m / 4'6"
Colour	Black / Red

*Range on full charge and flat ground

**Carrying capacity decreases to 113 kg / 18 stone with the maximum safe slope at 10°

Rascal Vantage X Profile:

The Rascal Vantage X is a compact, value for money 6 mph road scooter with extremely smooth user-friendly controls. With a comfortable Captain's seat as standard, infinitely adjustable delta-bar tiller steering and showerproof switches this highly manoeuvrable scooter features LED lighting and can adequately accommodate users up to 23.5 stone.

Local Rascal Mobility Retailer:

Main Features:

- Fold down backrest
- Flip-up width adjustable armrests
- Detachable storage basket
- Infinite adjustable delta tiller position
- Modern delta handlebars
- Emergency braking system
- Showerproof panel

**OPTIONAL
ACCESSORY**



Optional Extras:

- Standard all weather canopy - upgrade to a polycarbonate front shield for an additional cost
- Cane and crutch holder bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E - v10 Nov 2021

 electricmobility.co.uk