

Rascal 850

Performance scooters



Your new route to independence

RASCAL

Rascal 850

850 Specification:

| | |
|------------------------------|-----------------------|
| Type Class | C / 3 |
| Overall length | 127cm / 50" |
| Overall width | 65cm / 26" |
| Floor to seat height | 62cm |
| Seat width | 45.7cm / 18" |
| Seat depth | 43.1cm / 17" |
| Maximum carrying capacity | 136kgs / 21 stone |
| Battery voltage and capacity | 12V / 50Ah x 2 |
| Maximum speed | 12.9km/h / 8mph |
| Ground clearance | 15cm / 6" |
| Maximum kerb height | 10cm / 4" |
| Range* | up to 38km / 24 miles |
| The maximum safe slope | 9° |
| Turn-around width | 1.70m / 5'6" |
| Colour | Black / Red |

*Range of full charge and flat ground.
See full Tech Page on Website & Owner's Manual for details

Rascal 850 Profile:

Whether driving at 4mph on pavements or 8mph on roads, the Rascal 850 brings together compact design with Class 3 road legal specification. Available in black and with 50Ah batteries as standard you can travel in confidence and style. A fully adjustable sculpted sliding seat with added lumbar support combined with full suspension ensures a smooth ride, whatever the terrain.

Local Rascal Mobility Retailer:

Main Features:

- Five-position flip-up tiller
- Flip-up width adjustable armrests
- All round coil suspension
- Height adjustable headrest
- Backrest angle alteration
- Sturdy front bumper
- One wing mirror as standard

**OPTIONAL
ACCESSORY**

Optional Extras:

- Rascal armrest bag
- Cane & Crutch holder bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E.- v10 August 2025

 electricmobility.co.uk