


Vision

Performance scooters



Features **MaxGrip** Limited Slip Differential Gearbox



Your new route to independence

RASCAL

Vision

Vision Specification:

Type Class	C / 3
Overall length	154cm / 60.6"
Overall width	72cm / 28.3"
Floor to seat height	67-74.5cm (min-max)
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Maximum carrying capacity ¹	204kgs / 32 stone
Battery voltage and capacity*	12V / 75Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	12cm / 4.72"
Maximum kerb height	10cm / 4"
Range*- standard battery pack (75Ah)	up to 45km / 27.9 miles
Range*- upgraded battery pack (100Ah)	up to 60km / 37 miles
The maximum safe slope	10°
Turn-around width	2.21m / 7'3"
Colour	Black / Silver

*Range on full charge and flat ground

¹Max carrying capacity for 12V / 100Ah x2 battery pack - Ask your Retailer for more information

See full Tech Page on Website & Owner's Manual for details

Rascal Vision Profile:

The Rascal Vision offers a world of difference. More than you might expect from a Class 3 scooter, but everything you need to set the pace making life easier with just the right touch of style. The Vision combines a great modern look with uncompromising intelligent design and great suspension, **featuring MaxGrip Limited Slip Differential Gearbox**. Form and function come together to create something very special, with a well-padded seat providing good lumbar support and a charging point on its tiller, your days for adventure couldn't look better.

Local Rascal Mobility Retailer:

Main Features:

- Front coil and rear mono-shock suspension for a superior ride
- Gas strut tiller with easy access battery charging port
- Modern delta handlebars
- Modern LCD panel and twin mirrors
- Safety turning sensor, slows down scooter when turning at full speed
- Auxiliary USB power outlet for charging mobile phones etc
- Emergency braking system
- Modern LED lighting system

OPTIONAL
ACCESSORY

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- All weather canopy with a polycarbonate front shield



Max Grip Limited Slip Differential Gearbox Features:

- Improved grip, handling & safety
- Improved energy efficiency