Vision

Performance scooters



Features MaxGrip Limited Slip Differential Gearbox



Your new route to independence

Vision

Vision Specification:

Type Class	C/3
Overall length	154cm / 60.6"
Overall width	72cm / 28.3"
Floor to seat height	67-74.5cm (min-max)
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Maximum carrying capacity *	204kgs/32 stone
Battery voltage and capacity*	12V / 75Ah x 2
Maximum speed	12.9km/h/8mph
Ground clearance	12cm / 4.72"
Maximum kerb height	10cm / 4"
Range**- standard battery pack (75Ah)	up to $45 km / 27.9$ miles
	601 /07
Range**- upgraded battery pack (100Ah)	up to 60km / 37 miles
Range**- upgraded battery pack (100Ah) The maximum safe slope	up to 60km/3/ miles 10°
_	
The maximum safe slope	10°

^{**}Range on full charge and flat ground

Main Features

- Front coil and rear mono-shock suspension for a superior ride
- Gas strut tiller with easy access battery charging port
- Modern delta handlebars
- Modern LCD panel and twin mirrors
- Safety turning sensor, slows down scooter when turning at full speed
- Auxiliary USB power outlet for charging mobile phones etc
- Emergency braking system
- Modern LED lighting system

OPTIONAL ACCESSORY

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- All weather canopy with a polycarbonate front shield

Rascal Vision Profile

The Rascal Vision offers a world of difference. More than you might expect from a Class 3 scooter, but everything you need to set the pace making life easier with just the right touch of style. The Vision combines a great modern look with uncompromising intelligent design and great suspension, featuring MaxGrip Limited Slip Differential Gearbox. Form and function come together to create something very special, with a well-padded seat providing good lumbar support and a charging point on its tiller, your days for adventure couldn't look better.

Local Rascal Mobility Retailer:

- Max Grip Limited Slip Differential Gearbox Features:
- Improved grip, handling & safety
- Improved energy efficiency







[†] Max carrying capacity - batt pk (100Ah) Ask your Retailer *Upgrade to 12V / 100Ah x 2