

Frontier

Performance scooters



Your new route to independence



Frontier

Frontier Specification:

Type Class	C / 3
Overall length	143cm / 56"
Overall width	65cm / 26"
Floor to seat height	63 - 68cm
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity**	12V / 33Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range* - standard battery pack(33Ah)	up to 24km / 15 miles
Range*- upgraded battery pack (55Ah)	up to 40km / 25 miles
The maximum safe slope***	7°
Turn-around width	1.70m / 5'6"
Colour	Black / Red

*Range on full charge and flat ground **Upgrade to 12V / 55Ah x 2

***Maximum safe slope increases to 10° with a maximum carrying capacity of 106 kg / 16 stone with battery upgrade.

Rascal Frontier Profile:

The Rascal Frontier is an entry level 8 mph scooter, a proven performer with excellent legroom, comfort and dependability. It features leaf spring suspension, a delta-bar tiller with a comfortable swivel seat as well as showerproof switches and a large basket making the Frontier ideal for life on the move.

Available with an optional battery upgrade to 55Ah, that can give up to a 25 mile range.

Local Rascal Mobility Retailer:

Main Features:

- Flip-up adjustable armrests
- Leaf spring suspension
- Removable front basket
- Gas strut tiller adjustment with useful battery charging port
- Modern delta handlebars
- Showerproof panel and twin mirrors
- Large shock absorbing tyres for maximum comfort on the move

OPTIONAL ACCESSORY

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- Standard all weather canopy - upgrade to a polycarbonate front shield for an additional cost



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E. - v 9 Nov 2021

 electricmobility.co.uk