

# Ryley SL

Mid-wheel powerchairs



Your new route to independence



RASCAL

# Ryley SL

## Ryley SL Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	61cm / 24"
Floor to seat height Lowest to highest elevation position	56-75cm (min-max)
Seat width	47cm / 18.5"
Seat depth	48cm / 18.9"
Maximum carrying capacity	136kgs / 21.4 stone
Battery voltage and capacity	12V / 35Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 26.3km / 16.3 miles
The maximum safe slope	6°
Turn-around width	1.18m / 3'87"
Colour	Orange Sunset / Green Lightning

Single Heaviest component 65kg / 143.3lbs

\*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

## Rascal Ryley SL Profile:

The Rascal Ryley SL is a medium-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/ outdoor usage, where space is restricted & with durable full suspension. The Ryley SL has a range of up to 16.3 miles and delivers great performance with 35Ah batteries. This model includes a comfortable Captain's seat with a 7.5" column lift which delivers a smooth lift action and reaches a maximum floor to seat height of 75 cm when fully elevated.

Local Rascal Mobility Retailer:

## Main Features:

- 7.5 inch column lift
- Full suspension
- Outstanding manoeuvrability
- Flip-up fully adjustable armrests
- Height adjustable Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Dual density, puncture resistant drive wheels for maximum comfort on the move
- Maximum floor to seat height when fully elevated - 75 cm / 29.5"

### OPTIONAL UPGRADE

## Optional Upgrade:

- Swing-away legrests
- Elevating legrests
- Attendant Control
- Dual Control

### OPTIONAL ACCESSORY

## Optional Extras:

- Lighting Kit
- Armrest bag
- Rear Saddle bag
- Cane & Crutch holder bag
- Oxygen bottle bag



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E. - v 1Final-January 2023

 [electricmobility.co.uk](http://electricmobility.co.uk)