

Rivco SL



Your new route to independence

RASCAL 

Rivco SL

Rivco SL Specification:

Type Class	B / 2
Overall length	110cm / 43.3"
Overall width	70cm / 27.5"
Floor to seat height Lowest to highest elevation position	62-87cm (min-max)
Seat width	18.5", 20.5", 22"
Seat depth	18.9", 19.5", 19.5"
Maximum carrying capacity	182kgs / 28.6 stone
Battery voltage & capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 30.2 km / 18.8 miles
The maximum safe slope	6°
Turn-around width	1.6m / 5.2'
Colour	Blue Moon / Sunset Orange

*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

Rascal Rivco SL Profile:

The Rascal Rivco SL is a heavy-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, with durable full suspension & outstanding manoeuvrability. The Rivco SL has a range of up to 18.8 miles and delivers great performance with 55Ah batteries. This model includes a comfortable Captain's seat with a 10" scissor lift mechanism. The scissor lift delivers smooth and quiet lift action and reaches a maximum floor to seat height of 87 cm when fully elevated.

Local Rascal Mobility Retailer:

Main Features:

- 10 inch scissor lift mechanism
- Flip-up fully adjustable armrests
- Full suspension
- Height adjustable Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Maximum floor to seat height when fully elevated - 87 cm / 34"

OPTIONAL UPGRADE

Optional Upgrade:

- Swing-away legrests
- Elevating legrests
- Attendant Control
- Dual Control

OPTIONAL ACCESSORY

Optional Extras:

- Lighting Kit
- Armrest bag
- Rear Saddle bag
- Cane & Crutch holder bag
- Oxygen bottle bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v 1 Final-January 2023