

Rascal Razoo

Compact powerchairs



Agility made simple

Your new route to independence

RASCAL 

Razoo

Razoo Specification:

Type Class	B / 2
Overall length	89.6cm / 35.3"
Overall width	58.2cm / 22.9"
Floor to seat height	51.5 - 53 cm (min-max)
Seat width	47cm / 18.5"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21.4 stone
Battery voltage & capacity	12V / 22Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	3.5cm / 1.4"
Maximum kerb height	4cm / 1.6"
Range*	up to 16.8km / 10.5 miles
The maximum safe slope	6°
Turn-around width	97cm / 3'2"
Colour	White
Single Heaviest component	27.5kg / 60.6lbs

*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

Rascal Razoo Profile:

The Rascal Razoo is a slim & sleek everyday companion to assist and enhance your independence. Simplicity & functionality combine to create this delicate yet durable model with minimalism design a core component. The compact chassis makes this one of the most petite powerchairs on the market, with a quick release battery pack & easy to dismantle components. Transportation, storage and travel have never been easier!

This Rear-Wheel Drive powerchair offers superb traction & outstanding manoeuvrability with a desirable turning circle of just 97cm.

Local Rascal Mobility Retailer:

Main Features:

- Car transportable
- Flip-up fully adjustable armrests
- Swing-away joystick
- Transportable battery box
- Convenient joystick charging port
- Solid tyres all round
- Comfortable, rotating Captain's seat
- Flip-up footplate

OPTIONAL ACCESSORY

Optional Extras:

- Armrest bag
- Oxygen bottle bag
- Cane & Crutch holder bag
- Swing-away or elevating legrests also available



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E. - v 2 Final-Mar 2023

 electricmobility.co.uk