

ELECTRIC
MOBILITY®



The Tilmore

 **BRITISH
MADE**
Quality Craftsmanship

Independent Dual-Motor



Waterfall Back



Featuring Split Ottoman Extension




CosiChair



The Tilmore

Waterfall Back

The Tilmore Independent Dual-Motor Riser Recliner is designed to offer exceptional comfort from head to toe. Featuring independent control, from a seated position as the legs raise, the split ottoman footrest opens up & extends to give greater & longer hip to heel support. Adjustable, stitched Waterfall Back design offers the option to reposition the lumbar & head support maximising the personal comfort level for the head, neck & lower back. The handset has a 'Home Button' function which returns the chair to the standard seated position. The handset includes a practical USB port to charge a mobile phone or tablet. This Riser Recliner boasts both twin magazine & useful left & right hand control pockets.



Size Guide	Small	Medium	Large
Seat height	47cm / 18½"	50.8cm / 20"	53.3cm / 21"
Seat width	47cm / 18½"	50.8cm / 20"	54.6cm / 21½"
Seat depth	48.3cm / 19"	50.8cm / 20"	53.3cm / 21"
Arm height	16.5cm / 6½"	17.8cm / 7"	19.1cm / 7½"
Back height	71.1cm / 28"	73.7cm / 29"	76.2cm / 30"
Hip to Heel - raised footrest	94cm / 37"	97.8cm / 38½"	101.6cm / 40"
Overall width	81.3cm / 32"	83cm / 33"	88.9cm / 35"
Overall height	104cm / 41"	108cm / 42½"	111.8cm / 44"
Weight limit	25 stone / 158.8kg	25 stone / 158.8kg	25 stone / 158.8kg

Main Features & Benefits:

- Featuring Split Ottoman Extension
- Independent back rest & leg support
- Upto 25 stone weight limit
- 'Reset' Home Button
- Soft touch handset complete with integral USB port
- Antimacassar & arm covers included
- Tailored Lateral & Button Back available Made to Order
- 5 Year Guarantee on frames, mechanisms & all electrics
- Quality Craftsmanship, BRITISH MADE 



Como Asteroid



Como Bloom



Como Fawn



Como Spruce



Exclusively available from your local Cosi Chair Specialist

Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice. Images are for illustrative purposes only. Details are correct at time of going to print.

E. and O.E- v1 January 2024