

Your new route to independence



Vippi

Vippi Specification:

Type Class	B / 2
Overall length	100cm / 39.4"
Overall width	49.5cm / 19.5"
Floor to seat height	34-39cm (min-max)
Seat width	45cm / 17.7"
Seat depth	35cm / 13.7"
Maximum carrying capacity	130kgs / 20.4 stone
Battery voltage & capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	3.5cm / 1.3"
Maximum kerb height	5cm / 2"
Range*	up to 22.5km / 14 miles
The maximum safe slope	6°
Turn-around width	1.35m / 4'4"
Colour	Blue Moon / Green Lightning / Red Sky

Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- · Removable front basket
- Tiller adjustment
- · Solid tyres all-round
- Swivel seat



Optional:

Rascal arm rest bag



Single Heaviest component 15.4kgs / 33.9lbs *Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

Rascal Vippi Profile:

Trendy, transportable & fuss free, the Vippi is lightweight, versatile & compact. Designed & engineered to make every day easy & spontaneous, whether you are off on an adventure or running an essential errand. White stitching on the soft padded seat provides a premium finish complemented by a stylish black tiller with a fantastic pop of colour. Re-discover your power with this super scooter.

It dismantles in a flash with an easy take-apart mechanism deconstructing into 5 parts for storage & transportation with the single heaviest part weighing only 15.4kgs.

Local Rascal Mobility Retailer:

FCTRI



Specifications subject to change without prior notice Vippi LiFe Model with Lithium Ferro-Phosphate Battery available