

Vippi

Transportable scooters



Your new route to independence

RASCAL

Vippi

Vippi Specification:

Type Class	B / 2
Overall length	100cm / 39.4"
Overall width	49.5cm / 19.5"
Floor to seat height	34-39cm (min-max)
Seat width	45cm / 17.7"
Seat depth	35cm / 13.7"
Maximum carrying capacity	130kgs / 20.4 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	3.5cm / 1.3"
Maximum kerb height	5cm / 2"
Range* (Max user weight)	up to 15.6km / 9.7 miles
Range* (90kg user weight)	up to 18.5km / 11.5 miles
The maximum safe slope	6°
Turn-around width	1.35m / 4'4"
Colour	Blue Moon / Green Lightning / Red Sky
Single Heaviest component	15.4kgs / 33.9lbs

*Range on full charge and flat ground

*See full Tech Page on Website & Owner's Manual for details

Rascal Vippi Profile:

Trendy, transportable and fuss free, the Vippi is lightweight, versatile and compact.

Designed and engineered to make every day easy and spontaneous, whether you are off on an adventure or running an essential errand. White stitching on the soft padded seat provides a premium finish complemented by a stylish black tiller with a fantastic pop of colour. Re-discover your power with this super scooter.

It dismantles in a flash with an easy take-apart mechanism deconstructing into 5 parts for storage and transportation with the single heaviest part weighing only 15.4kgs.

Local Rascal Mobility Retailer:

Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Swivel seat



Dismantles easily into 5 parts



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v4 May 2026