

Veo Sport LiFe





† T&Cs apply

Veo Sport Life Specification:

Type Class	B/2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	39-44cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	130kgs / 20.5 stone
Battery voltage and capacity	25.6V / 6.2Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.4"
Maximum kerb height	5cm / 2"
Range*	up to 26.2km / 16.3miles
The maximum safe slope	6°
Turn-around width	1.40m / 4'6"
Colour	Neon Blue / Mulberry
Single Heaviest component	20.2kgs / 44.5lbs

Main Features:

- · British Made / Medical /Military grade **Lithium Ferro-Phosphate Battery**
- Car transportable
- All-round suspension
- · Stylish sports seat
- · Tiller charging point
- · LED front light
- Front bumper
- Dual compound puncture proof tyres
- Comfortable delta tiller



12 4Ah battery option

12.4Ah weighs just 4.7kg:

Relax & enjoy the ride in the knowledge that the Lithium Ferro-Phosphate (LiFePO4) Battery pack is the safest, greenest with an increased efficiency & performance to stand up to everyday life.

See full Tech Page on Website & Owner's Manual for details Rascal Veo Sport Life Profile:

*Range on full charge and flat ground

Offering peace of mind with a LiFetime Guarantee on the battery, it's all about the Stateof-The-Art Lithium Ferro-Phosphate Battery weighing just 4.7kg on this striking Rascal Veo Sport LiFe. Loaded with impressive features & modern, sporty design with class leading all-round suspension. With an Airline Compliant 12.4Ah battery, offering a range of up to 16.3 miles, plus padded delta tiller boasting a convenient charging point and luxury seat, detailed with unique white accent stitching, the Veo Sport LiFe is a real pleasure to ride! The Veo Sport LiFe's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension plus LED front light, this is the travel scooter with the edge.

Air travel & CAA Lithium Battery Legislation



Airline Compliant 12.4Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel & can be taken on a plane.

(Subject to Airline Regulation & Restrictions)

Now with 25Ah Lithium Ferro-**Phosphate Battery option**

- Range* upto 38.8km / 24.1miles
- 25Ah battery weighs 8kg
- Mass (weight) including batteries 51.6kg / 113.8lbs
- **NOT Airline Compliant**



† LiFetime Battery Guarantee T&Cs apply







Specifications subject to change without prior notice Images are for illustrative purposes only E. and O.E- v3 March 2024

