

Vista & Vista DX

Mid-size scooters



Your new route to independence

RASCAL



Rascal Vista & Vista DX Profile:

Discover the Rascal Vista and Vista DX mid-size pavement scooters. Benefiting from the quick release feature it dismantles into a few easy-to-lift sections for transporting in your car. The compact and slim-line style offers plenty of legroom yet it is still extremely robust with generous turning manoeuvrability, making it ideal for everyday use.

Each available in two eye catching colours with appealing features including easy access tiller charging and delta handle bar plus the all-round suspension offers added comfort. When taking apart, the split battery pack is also a real advantage. It comes with front and rear LED lights maximising visibility, giving you both confidence and support wherever you go. Available with solid puncture proof tyres only.



Enjoy day trips without the worry!

The Vista DX's increased 50Ah battery size offers an impressive range of up to 34.5 miles and the luxury Captain's seat makes it an excellent solution for active people who require a multi-purpose scooter and want to stay out for longer with friends and family.



Main Features Vista:

- All-round suspension
- Easy take apart action - car transportable
- Dual compound solid puncture proof tyres
- Comfort seat with folding backrest, adjustable arms, integrated swivel and slide functions
- Delta tiller with useful charging port
- Convenient front basket
- Class leading weight capacity
- Off-Board Battery Charging capability

Travel in comfort and style with the all-round suspension and swivel seating



Main Features Vista DX:

- All of the above features **PLUS**
- Captain's seat with headrest, folding backrest, adjustable arms, swivel and slide functions
- Range of up to 34.5 miles with 50A Battery capacity



Fully dismantles into 6 easy to handle parts for storage and transportation



Vista & Vista DX

Specification: Vista

Vista DX

Type Class	B / 2	✓
Overall length	118cm / 46.5"	✓
Overall width	58.5cm / 23"	✓
Floor to seat height	57cm - 62cm	✓
Seat width	44.5cm / 17.5"	45.7cm / 18"
Seat depth	43.1cm / 17"	43.1cm / 17"
Maximum carrying capacity	150kgs / 23.6 stone	140kgs / 22 stone
Battery voltage and capacity	12V / 35Ah x 2	12V / 50Ah x 2
Maximum speed	6.4km/h / 4mph	✓
Ground clearance	7.5cm / 3"	✓
Maximum kerb height	6cm / 2.3"	✓
Range*	up to 38km / 24 miles	up to 55.5km / 34.5 miles
The maximum safe slope	6°	✓
Turn-around width	1.7m / 5'8"	✓
Colour	Oxford Blue / Sky Red	Oxford Blue / Sky Red
Single Heaviest component (front Chassis)	25kgs / 55.12 lbs	✓

*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

Solid puncture proof tyres



Local Rascal Mobility Retailer: